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B Biodex Advantage - [Dynamometer Operation]

File View Window Wizard Help

Test Exercise Biofeedback Biofeedback Panel Control System Status Set ROM Limits Options Wizard

Verification Performed By : [Redacted]

Before ROM has been set:

1. Attach long shoulder attachment to dynamometer
2. Remove any handgrip
3. Click "ROM SET" button to bring up the ROM setup dialog (further instructions will be in the ROM dialog)

After ROM has been set:

4. Press "Mode" button and select Isometric mode
5. Move attachment to full horizontal (90 degrees), then press the "Start" button
6. Insert calibration weight until detent pin engages (torque = 50.0 ftlbs 67.8NM)
7. Click "Verify Calibration" button when ready
8. When finished, remove calibration weight and press the "Exit" button to leave the calibration verification screen

Torque is only updated when the "Verify Calibration" button is pressed Torque : [] N - m Result : []

Verify Calibration Report History EXIT

Mode: Setup Select Joint/Pattern End Stop: Cushion (1)

ROM Limits - % ROM		Speed	Torque	Contraction	Pause
Current ROM	Away	Away	Away	Away	Away
173	[]	[]	[]	[]	[]
Total ROM	Toward	Toward	Toward	Toward	Toward
0	[]	[]	[]	[]	[]
ROM SET	SET	SET	SET	SET	SET

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B Biodex Advantage - [Patient Selection]

File Record View Window Help

1 Patient 2 Protocol 3 Set ROM 4 Report

Open Add Patient Edit Save Cancel Del Test Del Patient Close

Last Name: Bonganha First Name: Valeria Gender: ☐ Male ☐ Female Involved: ☐ Right ☐ Left ☐ Both ☐ None

Height: 1.67 Weight: 59.5 (kgs) Birthdate: (M/d/yyyy) Phone: ID#: 3046

Address: Admission Date: Release Date: Diagnosis:

Test/Exercise Information

Date: 7/13/2015 3:18:22 PM Referral:

Clinician:

Notes:

Protocol Description: Isokinetic Unilateral Joint: Knee Pattern: Extension/Flexion

Pain Scale

0 1 2 3 4 5 6 7 8 9 10

No Pain Mild Discomforting Distressing Horrible Excruciating

Activate Report Generation Window NUM

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Biodex Advantage - [Protocol Definition]

File Record View Window Help

1 Patient 2 Protocol 3 Set ROM 4 Report

Protocol Favorites Add Edit Save Cancel Delete Linked Close

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Biodex Advantage - [Protocol Definition]

File Record View Window Help

1 Patient
2 Protocol
3 **Set ROM**
4 Report
Curve

Set Dynamometer Range of Motion

LEFT

Side Selected
LEFT ☒ RIGHT ☐

Use Previous ROM
AUTO SET ROM

Define New ROM
CLEAR LIMITS

AWAY Limit
SET RESET

TOWARD Limit
SET RESET

Legends
■ Patient ROM (Previous)
■ Patient ROM (Verified)

Total ROM:

Knee - Extension/Flexion
Chair Settings View Setup

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

Biodex Advantage - [261_1_uni_isom_conecc_1_1]

File View Window Help

1 Patient 2 Protocol 3 Set ROM 4 Report Curve

Print Preview Setup Close

Report Generation

☐ Progress Reports ☒ RIGHT ☒ LEFT (Involved)

Choose Options:

- ☐ Window isokinetic data
- ☐ Torque vs. Position Graph
- ☒ Use Metric Units
- ☐ Print as Unilateral
- ☐ Monochrome

Choose Report:

- General Evaluation
- Comprehensive Evaluation
- Rehab Session Summary

Print All Reports

Report Wizard Print Preview Print

Summary of TEST information

Subject: Valeria Bonganha
ID: 3046
Protocol: Unilateral
Mode: Isometric
Joint: Knee
Pattern: Extension/Flexion
Contraction: AWAY
Involved: Left
Session: 9/30/2011 3:43:23 PM

SET	SIDE	POSITION
1	LEFT	60
2	LEFT	60

Print

Printer: PDFCreator

Name: PDFCreator

Status: CutePDF Writer

Type: HP Deskjet 6940 series

Where: HP Deskjet 6940 series (Copy 1)

Comment: eDoc Printer

Print range: ☒ All ☐ Pages from: 1 to: ☐ Selection

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

Work Simulation	Closed Kinetic Chain	Lift Simulation	Dual position back
			

Características

- 17" Flat Panel LCD Touch Screen Monitor
- Multi-Mode operation; Isokinetic, Isometric, Isotonic, Reactive Eccentric and Passive
- Concentric speed up to 500 deg/sec
- Eccentric speed up to 300 deg/sec
- Concentric torque up to 500 ft-lb (680 Nm)
- Eccentric torque up to 400 ft-lb (444 Nm)
- Passive speed as low as .25 deg/sec
- - Passive torque as low as .5 ft-lb
- - Isotonic torque as low as .5 ft-lb.

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- Biodex System 4-Pro (Biodex Medical Systems, USA).