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- Click on Retrieve Subject Information for subject information from previous tests

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Select Test Profile:

- Select Density and TGV Models.
Click on Next > to continue

Density Model	Siri
Thoracic Gas Volume Model	Measured

Create Density Model

Subject Information More >>	
First Name	Sample
Middle Name	Test
Last Name	Subject
DOB	9/30/1976
Gender	Male
Height	190.0 cm

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
Test Complete:

- Ask subject to exit BOD POD
- Close BOD POD door
- Click on Next > to view test results

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Test Results:

- Modify entries using the provided entry boxes
- Click on Next > to save results

☐ Display / print results in lb and in

Test Results		More >>
% Fat	29.3 %	
% Fat Free Mass	70.7 %	
Fat Mass	15.923 kg	
Fat Free Mass	38.398 kg	
Body Mass	54.321 kg	
Estimated RMR	1055 kcal/day	
Estimated TEE	1350 kcal/day	
Daily Activity Level	Sedentary	

Subject Information
Test Profile

First Name	Sample
Middle Name	Test
Last Name	Subject
DOB	9/30/1976
Gender	Male
Height	190.0 cm
ID_1	3000
ID_2	4444
Ethnicity	General Population
Operator	operator
Test Date	
Test No.	

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BOD POD® Body Composition Tracking System Analysis

SUBJECT INFORMATION

NAME	Maria da Silva
AGE	34
GENDER	Female
HEIGHT	167.0 cm
ID_1	
ID_2	
ETHNICITY	General Population
OPERATOR	labfef
TEST DATE	April 16, 2015
TEST NUMBER	374

TEST PROFILE

DENSITY MODEL	Siri
THORACIC GAS VOLUME MODEL	Measured

BODY COMPOSITION RESULT

% FAT	19.1	%
% FAT FREE MASS	80.9	%
FAT MASS	11.138	kg
FAT FREE MASS	47.044	kg
BODY MASS	58.182	kg
BODY VOLUME	55.143	L
BODY DENSITY	1.0551	kg/L
THORACIC GAS VOLUME	3.867	L

OPERATOR COMMENTS

Body Fat: A certain amount of fat is absolutely necessary for good health. Fat plays an important role in protecting internal organs, providing energy, and regulating hormones. The minimal amount of "essential fat" is approximately 3-5% for men, and 12-15% for women. If too much fat accumulates over time, health may be compromised (see table below).

Fat Free Mass: Fat free mass is everything except fat. It includes muscle, water, bone, and internal organs. Muscle is the "metabolic engine" of the body that burns calories (fat) and plays an important role in maintaining strength and energy. Healthy levels of fat-free mass contribute to physical fitness and may prevent conditions such as osteoporosis.

BOD POD Body Fat Rating Table*

*Applies to adults ages 18 and older. Based on information from the American College of Sports Medicine, the American Council on Exercise, Exercise Physiology (4th Ed.) by McArdle, Katch, and Katch, and various scientific and epidemiological studies.

BODY FAT RATING	FEMALE	EXPLANATION
☐ Risky (high body fat)	> 40%	Ask your health care professional about how to safely modify your body composition.
☐ Excess Fat	30.1 - 40%	Indicates an excess accumulation of fat over time.
☐ Moderately Lean	22.1 - 30%	Fat level is generally acceptable for good health.
☒ Lean	18.1 - 22%	Lower body fat levels than many people. This range is generally excellent for health and longevity.
☐ Ultra Lean	15 - 18%	Fat levels often found in elite athletes.
☐ Risky (low body fat)	< 15%	Ask your health care professional about how to safely modify your body composition.

ENERGY EXPENDITURE RESULTS

Est. Resting Metabolic Rate (RMR) kcal/day	*Est. Total Energy Expenditure (TEE) kcal/day	Daily Activity Level
1259	1561	Sedentary
	1914	Low Active
	2191	Active
	2606	Very Active
*Est. TEE = Est. RMR x Daily Activity Level		

(See RMR Info Sheet for additional info)

Applies to adults ages 18 and older. Based on information from the Institute of Medicine (2002), Dietary Reference Intakes For Energy, Carbohydrate, Fiber, Fat, Fatty Acids, Cholesterol, Protein, And Amino Acids, Part I, pp93-206. Washington, D.C., National Academy of Sciences.

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Taxa Metabólica de Repouso – *Resting Metabolic Rate (RMR)*

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Nelson KM, Weinsier RL, Long CL, Schiti Y. Prediction of resting energy expenditure from fat-free mass and fat mass. Am J Clin Nutr; 1992; 56:848-56.

Fields DA, Hull, HR, Cheline, A.J, Yao M, Higgins, PB. Child specific thoracic gas volume prediction equations for air displacement plethysmography. Obes Res. 2004; 12: 1797–1804.

Demerath EW, Guo SS1, Chumlea WC, Towne B, Roche AF, Siervogel RM. Comparison of percent body fat estimates using air displacement plethysmography and hydrodensitometry in adults and children. Int J Obes. 2002; 26, 389–397.